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RT8 PRODUCT MANUAL

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RAYTOR INSTRUMENTS

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Dissolution Device and Solution Professional Provider

锐拓仪器
RAYTOR INSTRUMENTS

药物溶出度仪与溶出解决方案专业提供商
Dissolution Device and Solution Professional Provider

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RT800 Automatic Sampling Transdermal Diffusion System

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RT812 Automatic Sampling Transdermal Diffusion System

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After Sales Service & Technical Support

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Raytor Instruments is located in Shenzhen, the capital of Intelligent Manufacturing in China. It is a national high-tech enterprise and a ISO9001 certified enterprise.

Raytor Instruments is committed to the research and development, production, sales and service of drug dissolution instruments and solutions, with a number of technical patents.

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Raytor Instruments dissolution products include: flow-through cell dissolution system, reciprocating cylinder dissolution system, ordinary dissolution automatic sampling system, diffusion cell automatic sampling system, dissolution apparatus mechanical verification toolbox, dissolution medium preparation instrument, etc., which fully cover all kinds of dissolution release test devices in Pharmacopoeia. At the same time, Raytor Instruments can provide users with the overall solution of dissolution test.

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2nd Floor, Building 11, Nangang Second Industrial Park, Songbai Road, Nanshan District, Shenzhen

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Contact Number

400-0755-234 / 0755-23200354 / 18129931418

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Raytor Instruments Co., Ltd

邮箱(E-mail): sales@raytor.cn

网站(Website): www.raytor.cn

锐拓仪器
RAYTOR INSTRUMENTS

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Perfect Product Line

RT7

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Flow-Through Cell Dissolution System

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RT3

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Reciprocating Cylinder Dissolution System

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
Medium Preparation Degassing System

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RT2

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Dissolution Mechanical Calibration Series

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The performance of the instrument meets the requirements of Chinese Pharmacopoeia, United States Pharmacopeia, European Pharmacopoeia and other regulations, and the computerized system meets the requirements of 21CFR Part 11.

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In response to GMP inspection, prosecutors should observe the details as a breakthrough, so that the cleaning details of the equipment are impeccable.

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Excellent computer management system, domestic exclusive cloud management system, let the laboratory director control the equipment.

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We are not only producing instruments, we can provide professional technical support related to dissolution.

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The software system of Raytor dissolution apparatus is simple to operate and easy to use. It can store at least 200 dissolution test methods, and has at least three levels of permission level and audit track functions, data backup and PDF format export functions. It fully meets the requirements of data integrity of 21CFR PART 11.China GMP.EU GMP and other regulations.

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
At least 3 kinds of permission level can be allocated, including: system administrator, lab supervisor, general operator.

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The detailed rules of each permission level can be customized by customers according to their own risk assessment.

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At least 100 login accounts and passwords can be set up, and the user operation permission can be controlled by different levels.

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All levels of audit tracking, including: landing records, experimental records, operating records, cleaning records, etc.

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Audit track records can be retrieved, exported and printed.

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
The stand-alone version is equipped with a database system, the system storage capacity can meet at least 10 years of record storage.

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During the experiment, the test parameters can be real-time display, accurate and consistent electronic data can be generated synchronously.

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Real-time display of the UV/Vis detection results (for online UV / Vis detection).

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After the experiment, the system will automatically generate the experiment report.

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
Electronic data can be retrieved, exported, backed up and printed.

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The system database can be backed up and restored.

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Electronic data are protected from being tampered with or deleted.

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
Computerized System Technical Support Services

- 溶出仪计算机化系统验证服务。
Validation service of dissolution apparatus computerized system.
- 电子数据云系统服务器架设服务。
Electronic data cloud system server setup service.
- 操作系统可以进行升级和根据用户需求进行定制。
The operating system can be upgraded and customized according to user requirements.

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Raytor RT8 Transdermal Diffusion System meets the requirements of USP < 1724 >, which is used for drug transdermal release test of ointment, hard paste, coating, lotion, membrane, aerosol and other preparations, and studies the efficiency of drug transdermal absorption through the skin through the model of percutaneous absorption.

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Sampling error $< \pm 0.1$ mL. The function of rinse before sampling and cleaning after sampling, ensures sampling accuracy and avoids contamination or residue.

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Automatic Positioning Sampling Needle

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The sampling needle can be automatically positioned and tilted into the diffusion cell for automatic sampling. After that, it can collect the sample solution, and replenish the medium automatically.

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Medium Replenishment Cell (Constant Temperature)

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During the experiment, the medium to be replenished will be kept at a constant temperature. When replenishment is required, the sampling needle will automatically extract the medium from it and replenish it to the diffusion cell.

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8.4-inch touch screen, and pre installed with the Raytor operating system, meeting the data integrity requirements.

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Dry Heat Constant Temperature Stirring Device

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High precision temperature heating and stirring function. The automatic sampling experiment of 7 diffusion cells on both sides (14 in total) can be carried out simultaneously.

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Sample collection rack (7 × 12 position), suitable for liquid phase vials. During sampling, the sampling needle will automatically collect the sample solution into the liquid phase vial.

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Raytor improved Franz vertical diffusion cell is more portable and durable, and the unique sampling arm design makes the sampling process more convenient and smooth. The volume of the medium in the diffusion cell can be adjusted by using stirring bar with different heights.

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Users can set the method to control the system for automatic sampling, or manually use the pipette for sampling.

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
After sampling, the sampling needle will automatically move to the designated position of the sample collection rack to collect the sample solution into the liquid vial. It can support the sample collection of 7 diffusion cells on both sides at the same time.

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Each diffusion cell is equipped with an independent medium rewarming & replenishment cell, so that the temperature of the replenished medium can meet the requirements while effectively avoiding cross contamination between diffusion cell.

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It is used to measure the temperature of the medium in the vertical diffusion cell, and the temperature measurement error is less than $\pm 0.5^{\circ}\text{C}$.

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It is used to set the heating temperature and stirring speed, and control the start or stop of heating and stirring.

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More portable and durable, the unique sampling arm design makes the sampling process more convenient and smooth.

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The medium in the diffusion cell will be heated and mixed by dry heating and magnetic stirring.

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Install the vertical diffusion cell filled with medium on the dry heating stirring base. Start heating and stirring to preheat the medium. Use the temperature measuring probe to measure the temperature of the medium in each cell to ensure that the medium temperature has reached the specified value and remains stable.

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Adjust the heating temperature and stirring speed through the knob. The set value and current value of heating temperature and stirring speed will be displayed on the screen in real time.

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When dosing is required: Stop stirring, add the sample to the donor chamber assembly above the vertical diffusion cell, eliminate bubbles, restart stirring, and start the experiment. Conduct manual sampling through the sampling arm at the specified time point.

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Raytor provide vertical diffusion cells of different sizes, and the volume of media can be adjusted by using stirring bars of different heights.

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Raytor RT812 adopts automatic full sampling mode. During sampling, it will take out all the media in the diffusion cell and replenish fresh media of the same volume.

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Innovate "6+6" operation mode. The left and right diffusion cells (6 for each) can be set with different heating temperature and stirring speed, and two independent sets of percutaneous diffusion experiments can be conducted at the same time.

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Innovate the "6+6" sampling mode, which can simultaneously complete the sampling at different time points of two independent transdermal diffusion experiments.

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Preheating the medium (two kinds) for the experiment. Ensure that the medium entering the diffusion cell has been preheated to the temperature specified in the experiment.

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"6+6" operation mode. The left and right diffusion cells (6 for each) can be set with different heating temperature and stirring speed.

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
The volume of medium in the diffusion cell can be adjusted by using different size stirrer.

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Each cell is equipped with a temperature sensor.

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Manual sampling can be conducted through the sampling port. Or conduct full automatic sampling through the automatic sampling workstation.

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It is applicable to the transdermal diffusion test of photosensitive samples to avoid abnormal results caused by photodegradation.

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
"6+6" sampling mode, which can simultaneously complete the sampling at different time points of two independent transdermal diffusion experiments.

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Automatic full sampling. When sampling, take out all the media in the diffusion cell. After sampling, replenish fresh medium of the same volume.

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It can automatically clean the sampling pipeline and diffusion cell apparatus.

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The sample collection rack can be fitted with 1.5mL liquid phase vial and 10mL test tube. It can meet the sample collection requirements of various volumes.

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Thanks to the independently designed sample injection needle and sample collection rack. The automatic sampling workstation can realize the sampling of two groups of independent transdermal diffusion experiments at different time points and different sampling volumes.

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Automatic Sampling Transdermal Diffusion System

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Constant Temperature Water Bath & Stirring

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Preheat and keep the temperature of the media for the experiment. Preheating temperature and mixing speed can be set.

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Equipped with temperature sensor, which can monitor the temperature of heating in real time. The temperature sensor can be calibrated.

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In online mode, the RT812 system can be used to jointly control the medium thermostat apparatus, and the preheating temperature and stirring speed of the medium can be set in the transdermal diffusion test method.

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In offline mode, the preheating temperature and mixing speed can be set by adjusting the knob to meet the constant temperature mixing requirements of other experiments.

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

Raytor Instruments can provide comprehensive and professional validation and qualification services for your dissolution systems. All verification and validation services will be accompanied by complete and professional documents meeting regulatory requirements.

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The validation & qualification services we can provide include but are not limited to:

- 篮法、桨法 IQ/OQ/PQ
Basket Method, Paddle Method IQ/OQ/PQ
- 中国药典小杯法 IQ/OQ/PQ
ChP Mini Vessel Method IQ/OQ/PQ
- 桨碟法 IQ/OQ
Paddle over Disk Method IQ/OQ
- 转筒法 IQ/OQ
Dissolution Cylinder Method IQ/OQ
- 浸没池法 IQ/OQ
Immersion Cell Method IQ/OQ
- 流池法 IQ/OQ
Flow-Through Cell Method IQ/OQ
- 往复桶法、往复支架法 IQ/OQ
Reciprocating Cylinder and Holder IQ/OQ
- 扩散池法 IQ/OQ
Diffusion Cell Method IQ/OQ
- 计算机化系统验证
Computerized System Validation
- 溶出系统清洁验证
Dissolution System Cleaning Validation

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
Raytor Instruments after-sales service team can quickly respond to your instrument maintenance needs. Raytor Instruments provides services covering repair, preventive maintenance and requalification after repair, to keep your instrument at its best.

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Raytor Instruments can provide you with a full series of training on the operation of dissolution system. In addition, the professional technical team of Raytor Instruments can conduct deeper communication for specific technical details according to your experimental needs.

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Raytor Instrument has a strong application technology laboratory team, which can undertake the development of new dissolution methods from customers, such as flow-through cell method and reciprocating cylinder method. And assist the customer to investigate the abnormal dissolution results.

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
Electronic Data Cloud System Server

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The electronic data cloud system server is to connect all dissolution systems in multiple laboratories to realize the unified management and automatic backup of electronic data such as records, users and methods of all dissolution systems.

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
Raytor Instruments will share the latest dissolution research results in various national pharmaceutical technology seminars. You can pay attention to the participation arrangement of Raytor and communicate with our technical experts at the seminar site.

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Raytor Instruments will regularly hold a lecture hall online to introduce the latest dissolution cutting-edge technology. You can enter the live broadcast platform of the auditorium to watch the technical lectures and interact with technical experts online.

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The official account of WeChat public of Raytor Instruments will regularly share the application experiment cases, announce the arrangement of the lecture hall and technical seminar. You can pay attention to our official account for the latest information.

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The information, description and technical parameters in this document are subject to change without notice.

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